

Compare the Neopost IS-280 iMeter™ postage meter to the FP mymail™ MAX: Better speed. Better features. Better Total Value!



Costs comparison	Neopost IS-280	FP Mailing Solutions mymail™ MAX
"Mail Now, Pay Later" – postage advance	\$700 postage credit with no charge	Not available
Long term price protection	Yes	No
Postage funding account access	FREE with NeoFunds® Service ¹	No comparable service available
Discounts on USPS® First-Class rates	Yes	Yes
Discounts on USPS® commercial rates using NeoShip	Average savings of 11% ²	Not available
Proprietary battery needed for scale	None	Yes
Estimated ink cost per month	Same	Same
Postal rate change protection	FREE	FREE
Rental cost per month	\$19.75	\$15.95
Discounted shipping rates for packages, free address correction and tracking while meeting latest USPS® requirements	Available with NeoShip service for only \$3.00 per month	No comparable service available
Features comparison	Neopost IS-280	FP Mailing Solutions mymail™ MAX
Maximum parcel rates available	70 lbs	5 lbs
Screen display	2 3/4" x 1 5/8" (82% larger)	2 5/8" x 15/16"
Weighing platform dimensions	7 1/4" x 5 1/2" (almost 4 times larger)	3 1/2" x 3"
Thickness capacity for automated processing (Will my envelope fit?)	5/16"	1/4"
Rate shortcut keys and Rate Wizard	Yes	Not available
Letters per minute	Neopost IS-280 is up to 20% faster ³	
Letters per minute with weighing	Neopost IS-280 is over 100% faster ⁴	
Accounting departments	1 (10 optional)	3
USPS® eServices (delivery confirmation, return receipt, etc.)	Optional	Not available
Direct high speed (LAN) – easiest/fastest connection	Yes	Not available
PC connection	Yes	Yes
Envelope retention against misprinting (patented)	Yes	No
Ink cartridge prints until empty	Yes	No, stops after a certain time even if not empty

¹No interest if balance is paid in full each billing cycle

²Average savings on Commercial Based Priority Mail

³IS-280 processing speed of 18 letters per minute (lpm) versus mymail™ Max processing speed of 15 lpm

⁴Varied weight processing speed of 7 lpm for IS-280 versus 3 lpm for mymail™ MAX